

# Pros cons and challenges

Amplicons/marker genes/16S rRNA gene/microbiomics

Open & reproducible microbiome data analysis spring school  
Wageningen, The Netherlands, May 28-30, 2018



Turun yliopisto  
University of Turku



Gerben DA Hermes, PhD  
Laboratory of Microbiology,  
Wageningen University & Research

# Some learning goals

## Understand

- Critical steps in sample prep & data analysis
- Limitations (which conclusions can and can't you draw)
- Data (biological!) interpretation (day 2/3)



## Interpret literature sense & nonsense



Commentary


Suddenly everyone is a microbiota specialist!

S.A. Boers<sup>1</sup>, R. Jansen<sup>2</sup>, J.P. Hays<sup>1,\*</sup>

<sup>1</sup>) Department of Medical Microbiology and Infectious Diseases, Erasmus University Medical Centre Rotterdam, Rotterdam

<sup>2</sup>) Department of Molecular Biology, Regional Laboratory of Public Health Kennemerland, Haarlem, The Netherlands




# Literature & microbial biomarkers



Catherine A. Lozupone et al. Genome Res. 2013;23:1704-1714



# Microbial biomarkers for obesity




Walters et al. Meta-analyses of human gut microbes associated with obesity and IBD. FEBS Lett 2014


# Before sample prep/data generation

- How will you obtain fecal samples in your study?
- What do you want to do with it?
  - Functional (transcriptome -> quenching)
  - Metabolome
  - *In vitro* fermentations
- Freezing (how fast? What temperature? Freezer space)
- Pre-label tubes



# Typical amplicon sequencing workflow



# 1. DNA isolation



Kennedy NA, Walker AW, Berry SH, Duncan SH, Farquarson FM, et al. (2014) The Impact of Different DNA Extraction Kits and Laboratories upon the Assessment of Human Gut Microbiota Composition by 16S rRNA Gene Sequencing. PLOS ONE 9(2): e88982. <https://doi.org/10.1371/journal.pone.0088982>



Bacterial Family	Kit				Extraction Site		
	FastDNA 2 fold change	p	MoBio fold change	P	RINH fold change	p	Patients included
<i>Lachnospiraceae</i>	0.96 (0.74–1.25)	0.775	0.63 (0.49–0.81)	0.001	1.17 (0.95–1.44)	0.160	H3,H4,I1,I2
<i>Bacteroidaceae</i>	1.13 (0.79–1.63)	0.501	2.13 (1.49–3.05)	<0.001	1.09 (0.81–1.46)	0.561	H3,H4,I1,I2
<i>Ruminococcaceae</i>	0.94 (0.79–1.13)	0.524	1.32 (1.11–1.58)	0.005	0.95 (0.82–1.10)	0.516	H3,H4,I1
<i>Enterobacteriaceae</i>	1.08 (0.74–1.57)	0.695	0.61 (0.43–0.88)	0.016	0.85 (0.63–1.15)	0.311	I1,I2
<i>Sutterellaceae</i>	0.77 (0.18–3.37)	0.735	1.11 (0.26–4.69)	0.892	3.84 (1.18–12.46)	0.031	H3,H4,I1,I2
<i>Clostridiaceae</i>	1.00 (0.77–1.30)	0.976	0.46 (0.36–0.59)	<0.001	0.88 (0.71–1.08)	0.243	I1,I2
<i>Porphyromonadaceae</i>	1.46 (0.41–5.19)	0.560	4.03 (1.16–14.01)	0.035	0.70 (0.26–1.94)	0.502	H3,H4,I1,I2
<i>Erysipelotrichaceae</i>	1.21 (0.81–1.81)	0.361	0.32 (0.21–0.47)	<0.001	0.88 (0.64–1.22)	0.445	H3,H4,I1,I2
<i>Rikenellaceae</i>	0.35 (0.16–0.76)	0.016	0.72 (0.33–1.56)	0.418	0.65 (0.35–1.19)	0.181	H3,H4

RINH: Rowett Institute of Nutrition and Health.

Participants were excluded if all data points for that bacterial family were < 0.5%. Reference sample was from participant H3 using FastDNA method 1 and extracted at the Institute of Medical Sciences. Differences are shown as fold change with 95% confidence intervals.

doi:10.1371/journal.pone.0088982.t002



# Additional considerations

- ~~(Preferably)~~ **Don't** work with low DNA yield samples
  - Can “contaminate” other samples
- Use negative extraction controls
- Laboratory reagent specific microbiota

**Inherent bacterial DNA contamination of extraction and sequencing reagents may affect interpretation of microbiota in low bacterial biomass samples.**

[Glassing A](#)<sup>1</sup>, [Dowd SE](#)<sup>2</sup>, [Galandiuk S](#)<sup>3</sup>, [Davis B](#)<sup>4</sup>, [Chiodini RJ](#)<sup>5</sup>.

⊕ Author information



[this article](#) [search](#) [submit a manuscript](#) [register](#)

BMC Biol. 2014; 12: 87.  
Published online 2014 Nov 12. doi: [10.1186/s12915-014-0087-z](#) PMID: PMC4228153

**Reagent and laboratory contamination can critically impact sequence-based microbiome analyses**

Susannah J. Salter,<sup>1</sup> Michael J. Cox, Elena M. Turek, Szymon T. Galus, William O. Cookson, Minam F. Moffatt, Paul Turner, Julian Parkhill, Nicholas J. Loman, and Alan W. Walker<sup>2\*</sup>

[Author information](#) [Article notes](#) [Copyright and License information](#)



[this article](#) [search](#) [submit a manuscript](#) [register](#)

Microbiome. 2016; 4: 29.  
Published online 2016 Jun 23. doi: [10.1186/s40168-016-0172-3](#) PMID: PMC4917942

**Comparison of placenta samples with contamination controls does not provide evidence for a distinct placenta microbiota**

Abigail P. Lauder,<sup>1</sup> Aoife M. Roche,<sup>2</sup> Scott Sherrill-Mix, Aubrey Bailey, Alice L. Laughlin, Kyle Bittinger, Rita Leite, Michal A. Elovitz, Samuel Parry,<sup>3</sup> and Frederic D. Bushman<sup>1\*</sup>

[Author information](#) [Article notes](#) [Copyright and License information](#)

This article has been [cited by](#) other articles in PMC.

**Abstract**


Go to: 

# 2. Barcoded PCR:

## Primers & Region: Coverage, Resolution, bias


### 16S rRNA gene

- Target for detection & identification of bacteria.
- Ideal phylogenetic marker.
  - Universally distributed, functionally constant,
  - Sufficiently conserved, no horizontal transfer (?)
- More than 3,000,000 sequences in databases.
- Alternating variable and conserved sequence domains.



Nature Reviews | Microbiology

variability





base position



## 2. Barcoded PCR: Primers

### *In silico* Coverage

Check for your ecosystem of interest




- Americans didn't have *Bifidobacterium* up to ~2006
- No *Verrucomicrobia* in soil (actually 20%)

## 2. Barcoded PCR: Primers Resolution



Unequal sequence conservation

Be very careful with species level classification

Up to genus level identification is recommended!



# Be careful with functional interpretations!



"Luckily we found THIS bone, so we were able to reconstruct the whole creature..."

## 2. Barcoded PCR: Primers Bias



- Polymerase
- Primer-dimers
- Amplification bias: sequence of the bacteria themselves (high GC)
- Rate of chimera formation
- Sequencing specific sequence specific error rates
- Etc...

### Example

What percent of the product molecules contain an error after PCR (30 cycles) with different polymerases?

Polymerase	1 kb template	3 kb template
Phusion High-Fidelity DNA Polymerases (HF Buffer)	1.32%	3.96%
Phusion High-Fidelity DNA Polymerases (GC Buffer)	2.85%	8.55%
<i>Pyrrococcus furiosus</i> DNA polymerase	8.4%	25.2%
Taq DNA polymerase	68.4%	205.2%

The table above demonstrates the low error rate of Phusion DNA Polymerase. After 30 cycles of PCR amplifying a 3 kb template, only 3.96 % of the product DNA molecules contain 1 (nucleotide) error each. This means that 96.04 % of the product molecules are entirely error-free. In contrast, after the same PCR protocol performed with Taq DNA polymerase, every product molecule contains an average of 2 errors.




## 2. Barcoded PCR: Primers

### No universal primer (yet ?)

Table 1 | Comparison of sequencing technologies


	Read length	Maximum insert size	Run time (hours (h) or days (d))	Reads per run	Relative cost factor (per Mb)	Scale of reads per sample	Scale of samples per run	Raw error rate (%)			
								Total	Insertions	Deletions	Mismatches
ABI 3730	800 b	>1 Kb	2 h	96	100	$10^1$	$10^1$	0.001	<<0.1	<<0.1	<<0.1
454 FLX Titanium	300–400 b	800 b	9 h	$10^6$	1	$10^1$	$10^2$	1	<1	<<0.1	<<1
454 FLX+	500–600 b	1200 b	23 h	$10^6$	0.7	$10^1$	$10^2$				
Illumina GAIIx	76–101 b	500 b	6–9 d	$4 \times 10^8$	0.1	$10^3$ – $10^6$	$10^3$ – $10^4$	<1	<<1	<<1	<1
Illumina HiSeq 2000	101–151 b	500 b	9–15 d	$3 \times 10^9$	0.002	$10^3$ – $10^6$	$10^3$ – $10^4$				
Illumina MiSeq	36–151 b	500 b	4h–27 h	$10^7$	0.06	$10^4$	$10^2$				
PacBio	1100 b	>1 Kb	1.5 h	$3.5 \times 10^7$	1.5	$10^1$	$10^1$	15	13	1	1
IonTorrent	200 b	400 b	2–3 h	$1.5 \times 10^6$ – $3 \times 10^6$	0.4	$10^1$	$10^2$	2	1	1	<1

Kuczynski et al., Nature rev. 2011







## 2. Barcoded PCR: barcoding strategy




Amplification primers with annealing sites:



Amplification products:




Sequencing primers with annealing sites:





# Analysis pipeline



CLC Genomics Workbench

QIIME

Mothur


DADA2

Etc.....



# Analysis pipeline

## OTU picking



### Operational Taxonomic Unit

- traditionally clustered @ 97% -> 'species' proxy
- Now : 'sequences'
- Summarize to genus level

~88%

~92%

~95%



~97%

*Escherichia Coli* & *Salmonella* >98% similar

Horizontal gene transfer

# Analysis pipeline:

## OTU picking



97%


ered



**WAGENINGEN**  
UNIVERSITY & RESEARCH

# Analysis pipeline

## Taxonomic classification




Different methods of classification -> different name (scoring system)

- BLAST
- Kmer
- MisMatches

### 16s rRNA gene Databases

- Ribosomal Database Project (RDP-II)
- ~~Greengenes (13\_05)~~ (poor alignment/outdated)
- SILVA

# Challenges



Important!

Bias (might miss bacteria)  
Can overshadow biology

Important!

Bias (eg. coverage, amplification bias)  
Can overshadow biology

Important!

Bias (Sequencing error (type **and** machine specific))

Important!


Bias (OTU clustering)

Important!

Bias (classifier and database (different names))




The hard part

# Effects of all these different parameters in literature



Catherine A. Lozupone et al. *Genome Res.* 2013;23:1704-1714

# Microbial biomarkers



Walters et al. Meta-analyses of human gut microbes associated with obesity and IBD. FEBS Lett 2014



# Pros & cons

	Who is there?	Who is there, what can they do?	Who is there and who is doing what?		What are the enproducts?
	16S rRNA gene	Metagenomics	Metaproteomics	Meta transcriptomics	Metabolomics
What	Amplicons	DNA	RNA	Proteins	Metabolites
Cost/sample	€	€€€	€€€€	€€€	€€
Complexity	Medium	High	High	High	High
Taxonomic resolution	High	Very high	Medium	High	Low
Activity (precision)	-	-	Medium	High	High
Remarks	Bias associated with primers	Relatively low depth of analysis	Peptide spectrum matching is complex	RNA isolation challenging	No taxonomic information



# Keep in mind/challenges/interpretation

## OTUs

- Prokaryotic species level boundary is still hotly debated
- Species level concept for organisms with mobile genetic elements....
- Function & 16S are 'only' correlated (sometimes actually not that much)
- E coli: probiotic (*nissle*), pathogen (O157:H7 acute haemorrhagic colitis, Shiga-toxin) & everything in between (K12)
- Most bacteria are opportunists

# Recap